



GENESIS
STANDARD
STERILIZABLE
IN PLACE
SOLUTIONS



SOLARIS
BIOTECH SOLUTIONS

STANDARD STERILIZABLE IN PLACE SOLUTIONS

GENESIS

GENESIS is R&D Sterizable-In-Place Benchtop Fermenter/Bioreactor available from 7,5 up to 20 litres total volume. Automatic sterilization by steam or alternative through electric heaters (steam source not necessary) .

GENESIS is an ideal partner for microbial fermentation as well as animal, plant and insect cell cultivation. Typical applications includes the following:

- Education**
- Basic research**
- Scale-up and scale-down studies**
- Process development and optimization**

GENESIS can be used for:

- Biopharmaceutical**
- Biofuels research and manufacturing**
- Vaccines**
- Food and beverage biotechnologies**
- Bioremediation**
- Bioplastics**
- Cosmeceutical**
- Nutraceutical**

The best ratio
**Quality/
Capability/Price**
on the market

**WHY TO
INVEST**
IN THIS PRODUCT

**Automatic
sterilization**
through electrical heaters
(no need for an
external steam source)
or by steam



Benefits


Powerful/ Accurate **brushless motor**, from 1 to 2000 RPM.
Online absorbed Torques (Nm) and Power (W) measurements
obtaining an indirect density indication of the culture broth.

Sampling system




Illuminated side glass

Different gas mixing strategies with
up to 5 TMFC



External additional boxes parameters for
future PCS upgrade including dCO₂, Cell
Density, Weight, Peristaltic pumps, ect



SBC16: smart controller designed to provide
an high level of automated
management of the fermentation/
cultivation processes
Batch, Fed batch or continous processes

Modbus Digital
Hamilton sensors

Compact and modular PCS
(350x350x350 mm)

Double jacket (side-bottom)

Increased heat transfer
efficiency
It ensures optimal temperature
control and sterilization even at

N.4 assignable Watson Marlow pumps,
all speed controlled in entry level

Harvest valve in entry level
optionally SIP

Automatic sterilization by steam or
alternative through electric heaters




Modbus Hamilton sensors

Why a digital sensor?

Hamilton sensors (including Cell Density) has been integrated into Solaris PCS and Leonardo software giving the user the benefit of having a unique platform.

Fully compensated digital sensors, store and transmit all relevant sensor data, including calibration and diagnostic information directly to Solaris Leonardo software.



Sensor life
traceability

Reducing
background noise

pH

The electrolyte of the EasyFerm Bio sensors is prepressurized to prevent the diffusion of sample into the sensor. The Everef-F reference cartridge ensures that the reference electrolyte remains free of silver and precipitation of proteins.

dO2

The VisiFerm DO is the first optical oxygen sensor with integrated opto-electronics. The visiFerm requires less maintenance than a classical oxygen sensor as it does not have a mechanically sensitive membrane or a corrosive electrolyte.


ORP

The ORP sensor through a pre-pressurized reference electrolyte has a clog-free diaphragm.

The sensor ensures a stable measurement signals after steam sterilization, autoclavation and CIP cleanings with almost drift-free measurement.

Conductivity

All wetted conductivity sensor parts are FDA approved, can be cleaned easily and withstand CIP cleanings and autoclavations. The sensor shows a very good linearity over a broad measuring range.



ON LINE MEASUREMENT OF TOTAL CELL DENSITY




- Simple online measurement of cell growth
- Reliable values during the growth phase
- Early detection of process deviations

The Dencytee sensor performs online measurement of total cell density in solution. The sensor is based on optical density, which measures the turbidity of the cell suspension. The measurement is made at NIR (near-infra red) wavelengths so it is insensitive to changes in media color. All particles and molecules that scatter light at 880 nm will be detected, including living and dead cells as well as cell debris. This measurement is effective after inoculation when cells are expanding quickly but concentrations are low, making capacitance-based readings less reliable.

HOW IT WORKS

The Dencytee sensor emits light through a 5 mm window onto a light detector. Cells in suspension absorb and scatter light so less light is read by the detector. To compensate, the sensor increases the amount of light emitted by the light source to maintain a constant reading at the detector. By reading the amount of light that is increased at the light source, the Dencytee sensor can measure solutions with high cell densities.



ON LINE MEASUREMENT OF VIABLE CELL DENSITY




- Increase yield and lower production costs
- Detect changes in cell physiology with frequency scanning
- Precisely control harvesting for continuous culturing
- Early detection of process deviations

The Incyte sensor enables real-time, online measurement of viable cells in solution. The measurement is not influenced by changes in the media, microcarriers, dead cells or debris, and is designed for mammalian cell culture, yeast and high-density bacterial fermentation. Online measurement of viable cells makes it possible to detect events and respond in real time without sampling.

HOW IT WORKS

The Incyte measurement principle is based on capacitance. In an alternating electrical field, viable cells behave like small capacitors. The charge from these small capacitors is measured by the sensor and reported as permittivity (capacitance per area).



SALAS - Solaris Sterile Needle Free Additions System

SBC 16

NEEDLE
FREE

Genesis is supplied with **SALAS**, a 4 channels needle free additions system (Inoculums/ Feedings/pH corrective solutions/A.F solution).

EASY & QUICK
OPERATION

SALAS allows an easy and quick connection between the feeding solution and the vessel top lid.



USER-FRIENDLY SOFTWARE

The software is the user's best friend in experimental design planning and performing trial runs, as well as analyzing and optimizing media and parameters for cultivation.

The graphical user interface enables you to select the software functions intuitively. Data extracted are compatible with Windows Excel. However, Solaris has developed a platform where to easily and quickly manage fermentation data. This software is included in the fermenter supply and can be installed on unlimited number of client's PC or laptop.



Gas mixing

Various controller and hardware configurations enable aeration strategies using air, oxygen, nitrogen or a mixture of these to enrich the air. The mass-flow controller allows the exact flow rate control of individual gases. The flexible aeration options integrated in the bioreactor permit a wide range of different application giving to this system a substantial versatility.

- Thermal Mass Flow Controller in entry model
- Automatic gas mixing
- Gas mixing through TMFC and solenoid valves or numbers of TMFC
- Toro and sintered spargers



Data sheet

Vessel				
Solaris Code	Genesis 7.5	Genesis 10.0	Genesis 15.0	Genesis 20.0
Total Volume (liters)	7.5	10.0	15.0	20.0
Ratio D/H	1:2,5	1:2,5	1:2,5	1:2,5
Min. Working Volume (liters)	1.3	1.8	2.7	3.6
Max. Working Volume (liters)	5.6	7.5	11.25	15
Working temperature range	0-135°C			
Working pressure range	2 bar			
Design	Stainless Steel Jacketed Vessel			
Materials	Parts in contact with the culture AISI 316 L - other parts AISI 304			
Finishing	All parts in contact with the culture: Ra < 0,5 µm ; External: Ra < 0,6 µm Mirror polished			
Ports and Connections				
	Connection	Description		
Vessel lid	PG13	Antifoam		
	TC 3/4"	Safety valve		
	TC 3/4"	Gas-out		
	TK 3/4"	SALAS-Solaris Sterile liquid addition		
	TC 1"	Pressure probe		
Upper side wall	DN 52	Stirrer		
	TC 1/2"	Overlay gas inlet		
	TC 1/2"	Sparger		
	In gold	Sight glass		
	In gold	Sight glass		
Lower side wall	Hygenic socket	pH probe		
	Hygenic socket	dO probe		
	Hygenic socket	spare probe		
	Hygenic socket	spare probe		
	Temperature housing	PT100		
Vessel bottom	TC 3/4"	Harvest/sampling valve		
	TC 1/2"	Steam in		
	TC 1/2"	Water in		
	TC 1/2"	Jacket out		
Jacket in-out	1/2" G	Electric heaters		
	1/2" G	Electric heaters		
	1/2" G	Electric heaters		
Stirring				
Drive	Brushless Motor, Direct Assembly, 1-1500 rpm (bacterial), 1-500 (cell cultures)			
Power	208W (7.5-10L) ; 622W (15-20L)			
Impellers	Select from: Rushtons impellers , Marine Impellers, Pitched blade			
Thermoregulation				
Control	PID Control - Accuracy 0,1 °C			
	Jacket steam and electric heaters / cooling source			
Gas Control & Gas Mixing				
Sparger and overlay Gas Control	TMFC			
Gas Mixing (Air,CO ₂ ,O ₂ ,N ₂)	n.1 TMFC + n.4 solenoid valves, n° of TMFC			
Sparger type	Select from: Toro type (ring), syntered microbubbling both provided with 0,2 µm filter			
Exhaust	Condenser and 0,2 µm filter			
Controller				
Master Control Module	Dimensions Height: 350 mm Largeness: 350 mm Depth: 350 mm			
HMI with Leonardo software	23" touchscreen			

Controls

Temperature	
Sensor	PT100
Control system	Measuring resident in Leonardo 2.0 software
Control range	0 - 150°C
pH	
Sensor	Digital Hamilton sensor
Control system	Measuring resident in Leonardo 2.0 software
Control range	0 - 14
Operation temperature	0 - 130°C
Pressure range	0 - 6 bar
Actuator	Cascade to peristaltic pumps for the addition of acid/base solutions or gas (CO ₂)
dO ₂	
Sensor	Digital Optical Hamilton sensor
Control system	Measuring resident in Leonardo 2.0 software
Control range	0,05 - 300% air saturation
Operation temperature	-10 - 130°C
Pressure range	0 - 12 bar
Actuator	Cascade to RPM, Gas Control, feedings,ect
Antifoam/Level	
Sensor	Solaris sensor
Redox (ORP)	
Sensor	Digital Hamilton sensor
Control system	Measuring resident in Leonardo 2.0 software
Control range	±2000 mV
Operation temperature	- 10 -130°C
Pressure range	≤ 6 bar
Conductivity	
Sensor	Digital Hamilton sensor
Control system	Measuring resident in Leonardo 2.0 software
Control range	1 - 3000 µS/cm
Operation temperature	0 -130°C
Pressure range	0 - 20 bar
dCO ₂	
Sensor	Mettler Toledo sensor
Control system	Measuring resident in Leonardo 2.0 software
Control range	0,00-200% saturation
Operation temperature	-20.0-150°C
Pressure range	0 - 4 bar
Cell density	
Sensor	Hamilton-Fogale sensor
Control system	Measuring resident in Leonardo 2.0 software
Pressure range	0-3 bar (option 1) 0-10 bar (option 2)
Operation temperature	0-60°C (option 1) 0-80°C (option 2) (max. sterilization temperature 135°C)
Option 1	Total cell density based on turbidity (Two ranges: 10 ⁴ 5 to 10 ⁴ 8 mammalian cells/ml - 0.5 to 100 g/L dry weight)
Option 2	Viable cell density based on capacitance (Two ranges: 5x10 ⁴ 5 to 8x10 ⁴ 8 mammalian cells/ml - 5 to 200 g/L dry weight)
Weight	
Sensor	load cells
Control	Measuring resident in Leonardo 2.0 software
Peristaltic pumps	
WM 114	10-60 rpm
WM 313 FDM/D	45-350 rpm

Chiller

- Optionally GENESIS can be equipped with a chiller for heat removal from your culture minimizing lab water usage
- Using this system you don't need a water supply line in your lab
- Cost-effective cooling of fermenters
- Easy operation
- Refregerant level monitoring



Chiller data sheet

Working temperature range	-10°C / +40°C
Temperature stability	±0.5
Power consumption	0.7 kW
Filling volume range	2-8 L
Cooling output at 20°C measured with ethanol	0.25-0.60 kW
Cooling output at 10°C measured with ethanol	0.20-0.50 kW
Cooling output at 0°C measured with ethanol	0.15-0.36 kW
Cooling output at -10°C measured with ethanol	0.09-0.15 kW
Pump pressure max.	0.35-1.30 bar
Pump flow max.	16-35 L/min.



SOLARIS BIOTECHNOLOGY srl

Via Bachelet, 58 - 46047 Porto Mantovano
Mantova - Italy

Phone: +39 0376 408760

Fax: +39 0376 385108

Email: info@solarisbiotech.com

www.solarisbiotech.com